

ELASTICATED ADHESIVE BANDAGE 10 cm x 2.5 m non stretched



Product Code	34785
Unit of sale	box of 10 pcs
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	1848751
CND	M0304010301
EAN/UPC	7630016124976
PARAF	987474842

Description

Adhesive bandage size 2.5 m x 10 cm

(Length: 2.5 m non stretched, up to 4.5 m stretched)

Plaster bandages with hypoallergenic, synthetic adhesive containing zinc oxide, for permanent compression and functional dressings. Hypoallergenic, porous, air permeable with woven edges, can be sterilized by ethylene oxide. 100% cotton, latex free.

MADE IN SWITZERLAND.