

SINGLEFIX FINGER DRESSING - medium



Product Code	34801
Unit of sale	box of 100 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2450062
CND	M03030201
EAN/UPC	8033971553307

Description

SINGLEFIX® FINGER DRESSING - medium - box of 100

Finger dressing consisting of a cotton sheath and soft inner absorbent padding. Convenient wrist knot fixing. Easy to wear, it protects the injured finger, offering protection and comfort.

Multilingual instructions: GB, FR, IT, ES, PT, DE, HU, PL, BG, GR, RO, SE, SI, NO, NL, LT, LV, CZ, DK, FI, RU, TR, SK, EE, JP, Arabic.