

COTTON WOOL 250 g - Z folded



Product Code	34834
Unit of sale	box of 40 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2557596
CND	M010101
EAN/UPC	8023279348347
PARAF	986008528

Description

100% COTTON WOOL

- Cotton wool 250 g
- Type: Z folded
- Closing laces: Yes

Hydrophilic 100% cotton for sanitary and cosmetic use, hygiene, exudate absorption, non-injured skin disinfection.

Non-invasive Class I medical device.

Non sterile.

Multilingual pack: GB, FR, IT, DE, ES, FI, GR, HR, PL, PT, RO, SE.