

COTTON WOOL 1000 g - roll



Product Code	34838
Unit of sale	box of 10 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2557630
CND	M010101
EAN/UPC	8023279348385
PARAF	986008555

Description

100% COTTON WOOL

- Cotton wool 1000 g
- Type: roll
- Closing laces: No

Hydrophilic 100% cotton for sanitary and cosmetic use, hygiene, exudate absorption, non-injured skin disinfection.

Non-invasive Class I medical device.

Non sterile.

Multilingual pack: GB, FR, IT, DE, ES, FI, GR, HR, PL, PT, RO, SE.