


STERILE ADHESIVE DEVICE FOR CANNULA FIXATION



Product Code	34845
Unit of sale	box of 50 pcs
Minimum order	1
Type	Medical device
Class	I S

RDM (NSIS)	2514688
CND	M04010201
EAN/UPC	8015143004008
PARAF	987308350
GMDN	17787

Description

STERILE ADHESIVE DRESSING FOR CANNULA FIXATION - box of 50 pcs

Designed for wide applications.

Hypoallergenic acrylic adhesive, dermatologically tested.

Size 8 x 5.8 cm.

Wound pad size 2 x 1.5 cm.

Multilingual box: GB, FR, IT, ES, DE.

Medical Device Class I.

Made in Italy.