



Product Code	34866
Unit of sale	box of 50 pcs
Minimum order	1
Type	Medical device
Class	I S

RDM (NSIS)	2511047
CND	M04010102
EAN/UPC	8015143003742
PARAF	987682174

## Description

Sterile hypoallergenic adhesive dressing with central non-wound-sticking absorbant pad.  
Sterile impermeable protection of surgical wounds, superficial wounds, cuts and abrasions.

- BWT dressing hospital line
- Size: 5 x 7.2 cm
- Wound pad: 2.8 x 3.8 cm
- Box of 50 pcs

Multilingual box: GB, FR, IT, ES, DE.  
Made in Italy.

- Waterproof
- BWT material: transpiring, waterproof, thin
- Round corners avoid dressing from coming off easily
- Polyurethane film acts as an antibacterial barrier