

TEGADERM 3M IV TRANSPARENT FILM DRESSING

adhesive-free window 12x12 cm



Product Code	34891
Unit of sale	box of 25 pcs
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	2541810
CND	M04010202
GTIN	00707387770389

Description

3M™ TEGADERM™ FILM

- 3M Code: 1668
- Product: Tegaderm™ - 12x12 cm with adhesive-free window - sterile

Transparent, sterile waterproof dressing.
Highly breathable.
Frame style for easy application also available in non-sterile roll.