

TEGADERM PORT 3M CHG IV 12x12 cm



Product Code	34892
Unit of sale	box of 25 pcs
Minimum order	1
Type	Kit of medical devices
Class	

RDM (NSIS)	144360
GTIN	30707387802210

Description

3M™ TEGADERM™ CHG IV SECUREMENT DRESSING

3M code: 1665R
Product: Tegaderm™ CHG - 12x12 cm

Integrated I.V. dressing combines infection reduction with site visibility, catheter securement, and breathability. Proven to reduce CRBSI and vascular catheter colonization. Immediate and continuous antimicrobial protection through integrated chlorhexidine gluconate gel pad, without requiring moisture to activate it. Designed to minimize catheter movement, it can be worn up to 7 days. Waterproof and sterile barrier to external contaminants (liquids, bacteria, and viruses).