

# TEGADERM 3M + PAD - 9x10 cm - sterile



Product Code	34896
Unit of sale	box of 25 pcs
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	2393907
CND	M04010102
EAN/UPC	8711428086928
PARAF	977814831

## Description

3M medication products - Sterile transparent film dressing - 3586  
TEGADERM™ - PAD 9 x 10 cm - sterile  
Minimum order: 1 box of 25 pcs

Transparent, sterile waterproof dressing. Highly breathable.  
Frame style for easy application.  
With a non-adherent absorbent PAD.  
Ideal for surgical wounds.