

TEGADERM 3M + PAD - 9x15 cm - sterile



Product Code	34897
Unit of sale	box of 25 pcs
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	2393909
CND	M04010102
EAN/UPC	8711428086904
PARAF	977814843

Description

3M medication products - Sterile transparent film dressing - 3589
TEGADERM™ - PAD 9 x 15 cm - sterile
Minimum order: 1 box of 25 pcs

Transparent, sterile waterproof dressing. Highly breathable.
Frame style for easy application.
With a non-adherent absorbent PAD.
Ideal for surgical wounds.