

TEGADERM 3M - 6x7 cm - sterile



Product Code	34900
Unit of sale	box of 100 pcs
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	2542639
CND	M04040601
EAN/UPC	8711428065510
PARAF	977814870

Description

3M™ TEGADERM™ FILM

- 3M Code: 1624W
- Product: Tegaderm™ - 6x7 cm - sterile
- Minimum Order: Box of 100

Transparent, sterile waterproof dressing.
 Highly breathable.
 Frame style for easy application also available in non-sterile roll.