



Product Code	34901
Unit of sale	box of 50 pcs
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	2542641
CND	M04040601
EAN/UPC	4054596115422
PARAF	977814882

Description

3M™ TEGADERM™ FILM

- 3M Code: 1626W
- Product: Tegaderm™ - 10x12 cm - sterile
- Minimum Order: Box of 50

Transparent, sterile waterproof dressing.
Highly breathable.
Frame style for easy application also available in non-sterile roll.