

TEGADERM 3M - 15x20 cm - sterile



Product Code	34903
Unit of sale	box of 10 pcs
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	2542650
CND	M04040601
EAN/UPC	8711428065411
PARAF	977814906

Description

3M™ TEGADERM™ FILM

- 3M Code: 1628
- Product: Tegaderm™ - 15x20 cm - sterile
- Minimum Order: Box of 10

Transparent, sterile waterproof dressing.
Highly breathable.
Frame style for easy application also available in non-sterile roll.