

STERI-STRIP 3M - 100 x 6 mm



Product Code	34911
Unit of sale	50 bags of 10 pcs
Minimum order	1
Type	Medical device
Class	I S

RDM (NSIS)	2183490
CND	H900102
PARAF	903666459
GTIN	30707387804665

Description

3M STERI-STRIP™ - SKIN CLOSURE STRIPS

REINFORCED STERI-STRIP™

With polyester reinforcing filaments for higher tensile strength, reduced infection and scarring risk and better cosmetic outcomes.

- STERI-STRIP™ - 6 x 100 mm
- Reinforced
- Sterile
- 3M Code: R1546
- Minimum order: 50 bags of 10 strips

Adhesive skin closures made of non-woven backing coated with a pressure-sensitive hypoallergenic adhesive. Sterile, breathable and comfortable to wear.