

STERI-STRIP 3M - 50 x 12 mm



Product Code	34913
Unit of sale	50 bags of 6 pcs
Minimum order	1
Type	Medical device
Class	I S

RDM (NSIS)	2183495
CND	H900102
EAN/UPC	8023279349139
PARAF	977814971
GTIN	30707387804719

Description

3M STERI-STRIP™ - SKIN CLOSURE STRIPS

REINFORCED STERI-STRIP™

With polyester reinforcing filaments for higher tensile strength, reduced infection and scarring risk and better cosmetic outcomes.

- STERI-STRIP™ - 12 x 50 mm
- Reinforced
- Sterile
- 3M Code: R1549
- Minimum order: 50 bags of 6 strips

Adhesive skin closures made of non-woven backing coated with a pressure-sensitive hypoallergenic adhesive. Sterile, breathable and comfortable to wear.