

STERI-STRIP 3M - 100 x 12 mm - Mini bag



|               |                |
|---------------|----------------|
| Product Code  | 34918          |
| Unit of sale  | bag of 6 pcs   |
| Minimum order | 1              |
| Type          | Medical device |
| Class         | I S            |

|            |               |
|------------|---------------|
| RDM (NSIS) | 2183491       |
| CND        | H900102       |
| EAN/UPC    | 8023279349184 |

Description

3M medication products - Sterile skin closure  
STERI-STRIP™ - 100 x 12 mm  
Minimum order: 1 bag of 6 strips

Sterile skin closure strips that are made of non-woven backing coated with a pressure-sensitive, hypoallergenic adhesive and reinforced with polyester filaments for added strength.