

TEGADERM 3M I.V. FIXING STRIPS 8.5x10.5 cm



Product Code	34946
Unit of sale	box of 50 pcs
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	2265716
CND	M04010202
EAN/UPC	4001895928512
PARAF	977815176

Description

3M TEGADERM™ IV - STERILE - 8.5 x 10.5 cm - for central lines
3M code: 1635

TEGADERM™ IV

A cost-effective solution for peripheral and central vascular lines with strong hold and gentle removal.

A complete range of sterile, waterproof, highly-breathable transparent film dressings to protect IV sites from water, virus and bacteria. With antifriction backing and reinforced non-woven borders for longer wearing time and minimized catheter dislodgement risk.