

TEGADERM 3M I.V. ADVANCED 6.5x7 cm for peripheral lines



Product Code	34947
Unit of sale	box of 100 pcs
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	2564113
CND	M04010202
EAN/UPC	4046719981120

Description

TEGADERM™ IV ADVANCED

A reinforced, 100% waterproof, advanced solution for better securement and longer wearing time.

- 3M code: 1683
- 6.5 x 7 cm for peripheral lines
- Sterile

A complete range of sterile, waterproof, highly-breathable transparent film dressings to protect IV sites from water, virus and bacteria. With antifriction backing and reinforced non-woven borders for longer wearing time and minimized catheter dislodgement risk.