

TEGADERM 3M FOAM 5x5 cm - non adhesive



Product Code	34953
Unit of sale	box of 10 pcs
Minimum order	1
Type	Medical device
Class	II B

RDM (NSIS)	2671221
CND	M04040601
GTIN	00707387807184

Description

3M™ TEGADERM™ HP FOAM DRESSINGS

3M code: 90600
Product: Tegaderm™ HP foam
Size: 5x5 cm
Minimum order: Box of 10
Type: Non-adhesive dressing

3M™ TEGADERM™ HIGH PERFORMANCE FOAM
Sterile, multi-layer foam dressing that combines high absorbency with high breathability to reduce the risk of wound maceration. Well suited for changing levels of exudate to maintain moisture balance for optimal wound healing. Highly breathable and waterproof, also available in heel/elbow shape with adhesive border and simplified application system.