

TEGADERM 3M FOAM 10x10 cm - non adhesive



Product Code	34954
Unit of sale	box of 10 pcs
Minimum order	1
Type	Medical device
Class	II B

RDM (NSIS)	2671223
CND	M04040601
EAN/UPC	0707387807207
PARAF	977815214

Description

3M™ TEGADERM™ HP FOAM DRESSINGS

3M code: 90601
Product: Tegaderm™ HP foam
Size: 10x10 cm
Minimum order: Box of 10
Type: Non-adhesive dressing

3M™ TEGADERM™ HIGH PERFORMANCE FOAM

Sterile, multi-layer foam dressing that combines high absorbency with high breathability to reduce the risk of wound maceration. Well suited for changing levels of exudate to maintain moisture balance for optimal wound healing. Highly breathable and waterproof, also available in heel/elbow shape with adhesive border and simplified application system.