

# TEGADERM 3M FOAM 14x14 cm - adhesive



Product Code	34957
Unit of sale	box of 5 pcs
Minimum order	1
Type	Medical device
Class	II B

RDM (NSIS)	2238110
CND	M04040602
EAN/UPC	4064035037826
PARAF	977815238

## Description

### 3M™ TEGADERM™ HP FOAM DRESSINGS

3M code: 90619

Product: Tegaderm™ HP foam adhesive

Size: 14x14 cm

Minimum order: Box of 5

Type: Adhesive dressing

Indication: For heel and elbow

### 3M™ TEGADERM™ HIGH PERFORMANCE FOAM

Sterile, multi-layer foam dressing that combines high absorbency with high breathability to reduce the risk of wound maceration.

Well suited for changing levels of exudate to maintain moisture balance for optimal wound healing.

Highly breathable and waterproof, also available in heel/elbow shape with adhesive border and simplified application system.