

TEGADERM 3M CHG 10x12 cm



Product Code	34960
Unit of sale	box of 25 pcs
Minimum order	1
Type	Medical device
Class	III

RDM (NSIS)	2395614
CND	M04010299
GTIN	30707387782185

Description

3M™ TEGADERM™ CHG IV SECUREMENT DRESSING

3M code: 1658R
Product: Tegaderm™ CHG - 10x12 cm
Minimum order: box of 25

Integrated I.V. dressing combines infection reduction with site visibility, catheter securement, and breathability. Proven to reduce CRBSI and vascular catheter colonization. Immediate and continuous antimicrobial protection through integrated chlorhexidine gluconate gel pad, without requiring moisture to activate it. Designed to minimize catheter movement, it can be worn up to 7 days. Waterproof and sterile barrier to external contaminants (liquids, bacteria, and viruses).