

CAVILON 3M HIGH PROTECION BARRIER FILM LIQUID

0.7 ml



Product Code	34978
Unit of sale	box of 4 pcs
Minimum order	1
Type	Medical device
Class	II A

UK-REP	Yes
--------	-----

RDM (NSIS)	2308041
CND	M04041503
PARAF	950833741

Description

Cavilon™ Advanced skin protectant - liquid 0.7 ml

3M™ code: 5051G4P

Minimum order: box of 4

3M™ Cavilon™ Advanced Skin Protectant

Highly durable, ultra-thin, transparent barrier which protects against the harmful effects of bodily fluids. Liquid solution with single-use applicator that reduces the potential for cross-contamination and IAD pain, creating a waterproof barrier film.

A complete range of products to protect skin from body fluids and friction.

Cavilon™ products create on the skin a highly breathable, long lasting transparent coating to prevent maceration and trauma.

Alcohol-free, hypoallergenic and non-irritating, can be used on broken skin and under any tapes and dressing.