

**CUROS STOPPER DISINFECTING CAPS - open female
Luer connector - 50 strips of 5**



Product Code	34982
Unit of sale	box of 250 pcs
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	2316282
CND	A0706
EAN/UPC	4054596459052

Description

3M™ CUROS™ DISINFECTING CAPS

- 3M code: CSV5-250R
- Curos™ Stopper™ caps
- Connector: open female Luer
- Sale formats: 50 strips of 5
- box of 250

3M™ CUROS™ DISINFECTING CAPS

3M™ Curos caps disinfect in 1 minute and protect needle connectors and Luer to help reducing the risk of contaminants from entering the catheter post-insertion, for up to 7 days if not removed.
Caps twist on easily and stay securely in place, preventing contained alcohol (70% isopropyl alcohol) from spreading within patient's circulatory system.
For single use only.
Suitable for needle-free connectors (34980), male Luer (34981) and open female Luer (34982).