

ADAPTIC DIGIT 3M NON ADHERENT DRESSING - for finger - large - sterile



Product Code	34995
Unit of sale	box of 10 pcs
Minimum order	1
Type	Medical device
Class	I S

UK-REP	Yes
--------	-----

RDM (NSIS)	2301032
CND	M02030201
GTIN	25051978001444

Description

- 3M™ ADAPTIC™ DRESSING - NON ADHERENT - sterile
- 3M Code: MAD023
 - ADAPTIC™ DIGIT - for finger - large
 - Minimum order: box of 10

Non-adherent dressing that combines a low tack silicone coated material with a tubular bandage specifically designed for the treatment of finger and toe wounds.

Conformable and comfortable dressing, which allows for flexibility and maximum movement of the fingers and toes.

Easy to apply, helps minimize trauma associated with dressing changes.

It can be left on for several days.