

DENTAL COTTON ROLLS - 10 pack of 1000



| | |
|---------------|------------------|
| Product Code | 35000 |
| Unit of sale | box of 10000 pcs |
| Minimum order | 1 |
| Type | Medical device |
| Class | I |

| | |
|--------|-----|
| UK-REP | No |
| CH-REP | Yes |

| | |
|------------|---------------|
| RDM (NSIS) | 2558968 |
| CND | Q019004 |
| EAN/UPC | 8023279350005 |
| PARAF | 987472533 |
| GMDN | 31815 |

Description

COTTON ROLLS 1 x 3.8 cm

100% cotton interior, soft, pliable, non-linting, non-irritating, highly absorbent up to 6 times their weight. Will not adhere to mucous membranes, maintains shape better when wet, non sterile.
Minimum order: carton of 10 boxes of 1,000 pcs.