

NON WOVEN GAUZE 36x40 cm - box of 12



Product Code	35027
Unit of sale	40 boxes
Minimum order	1
Type	Medical device
Class	I S

RDM (NSIS)	2294123
CND	M0202010101
EAN/UPC	3800232426176
GMDN	47107

Description

NON WOVEN FABRIC GAUZE WITHOUT X-RAY

- Size: 36x40 cm
- Minimum order: carton of 40 box of 12 pcs

Non invasive device to be used as a mechanical barrier, for compression or for absorption of exudates.
Intended to manage the micro environment of a wound.
Maximum absorbance, it does not adhere to the wound.
Made of viscose and polyester.
Medical device class risk I

Individually packed.
Made in Europe.