

COTTON GAUZE SWAB 5x5 cm - X-ray - sterile



Product Code	35036
Unit of sale	box of 350 pcs
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	2050629
CND	M0201020201
EAN/UPC	3800232420426

Description

COTTON GAUZE SWABS - sterile

Plies: 12

Threads: 17

Gauze swabs - sterile: 5x5 cm - X-ray

Minimum order: dispenser of 350 (70 packs of 5 swabs)

100% cotton gauze, 17 or 20 threads cm², according to European Pharmacopoeia and EN 14079.

Multilingual label and instructions: GB, FR, IT, ES, PT, DE, PL, SE, GR.