

COTTON GAUZE SWAB 7.5x7.5 cm - sterile



Product Code	35037
Unit of sale	box of 250 pcs
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	2050755
CND	M0201020101
EAN/UPC	3800232420433

Description

COTTON GAUZE SWABS - sterile

Plies: 12
Threads: 17
Gauze swabs - sterile: 7.5x7.5 cm
Minimum order: dispenser of 250 (50 packs of 5 swabs)

100% cotton gauze, 17 or 20 threads cm², according to European Pharmacopoeia and EN 14079.

Multilingual label and instructions: GB, FR, IT, ES, PT, DE, PL, SE, GR.