

COTTON GAUZE SWAB 10x20 cm - X-ray - sterile



Product Code	35042
Unit of sale	box of 100 pcs
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	2420745
CND	M0201020201
EAN/UPC	3800232427333
GMDN	15066

Description

COTTON GAUZE SWABS - sterile

Plies: 16
Threads: 20
Gauze swabs - sterile: 10x20 cm - X-ray
Minimum order: dispenser of 100 (10 packs of 10 swabs)

100% cotton gauze, 17 or 20 threads cm², according to European Pharmacopoeia and EN 14079.

Multilingual label and instructions: GB, FR, IT, ES, PT, DE, PL, SE, GR.