



Product Code	35058
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	Yes
RDM (NSIS)	2414667
CND	Z1203020408
EAN/UPC	6934440102311
UMDNS	17148
PARAF	987746866

## Description

### OXYFIT™ FINGER PULSE OXIMETER

Professional, soft and wearable fingertip oximeter for continuous monitoring of SpO2 and Pulse Rate during sleep or daily use. It measures, displays and stores adult results anytime, anywhere. Its unique duckbill probe has a soft chamber to check values in a comfortable and precise way.

- bright OLED display to show clear results
- built-in memory to store up to 4 sessions of record
- bar graph & pulse waveform display
- real time data can be recorded on free ViHealth App
- built-in lithium battery - automatic ON /OFF function
- sound alarms to report irregular readings.

Supplied with charging cable and user manual in GB, FR, IT , ES .

Size: 38x30x38 mm  
Weight: 28 g

## Technical Specifications

SpO2 range: 70% to 100%  
 SpO2 accuracy: 80-100%:±2%, 70-80%:±3%  
 Pulse Rate range: 30 to 250 bpm  
 Pulse Rate accuracy: ±2 bpm or ±2%, whichever is greater  
 Beep reminder source: low oxygen level, high/low pulse rate  
 Recorded parameters: oxygen level, pulse rate  
 Record intervals: 4s  
 Built-in memory: 4 records, up to 1 hours for each  
 Built-in battery: rechargeable lithium-polymer battery; 3.7Vdc  
 Charger: 5VDC, Max. 80mA  
 Charging time: 2-3 charging hours  
 Battery life: 12-14 hours continuous working  
 Bluetooth: Bluetooth 4.0 BLE  
 Protection: electric shock: Type BF electromagnetic: Group I, Class B dust and water resistance: IP22