



Product Code	35084
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II B

RDM (NSIS)	451807
CND	Z1203020408
EAN/UPC	0833166004847
GMDN	17148

Description

NONIN ONYX® VANTAGE 9590 PULSE OXIMETER - blue

The only finger pulse oximeters with scientifically proven accuracy in motion and low perfusion. Tested for use on fingers, thumbs, or toes.

- Accurate - Scientifically proven performance in the widest range of patient populations and settings.
- Durable - Protects against dropping and water spills; allows for thousands of uses in any environments and exceed IP32 water ingress testing.
- Efficient - Up to 6000 spot checks or 36 hours of continual operation on two AAA batteries.
- Safe - Built in the USA with lead-free and latex-free construction.
- Versatile - suitable for pediatric to larger adult patients. Wide range of digit thicknesses from 8 mm to 25.4 mm.
- Warranted - Industry-leading 4-year warranty.
- Cost Effective - Proven accuracy, quality and durability make the Onyx Vantage a smart purchase for the long term.

Multilingual CD manual: GB, FR, IT, DE, ES, PT, GR, NL, PL, NO, FI, DK, SE.

Large bright LED Display. SpO2 and pulse rate can be read from any angle, day or night.
Lead (PB) free and latex free

Technical Specifications

Pulse rate range: 18-321 bpm \pm 3 digit (20-250 bpm)
 SpO2 range: 0-100% \pm 2 digit (70-100%)
 Low perfusion SpO2: 70-100% \pm 2 digit
 Low perfusion pulse rate: 40-240 bpm \pm 3 digit