

OXY-4 FINGER OXIMETER - blue - GB,FR,IT,DE,ES,PT,GR,PL



Product Code	35091
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	1844279
CND	Z1203020408
EAN/UPC	8023279350913
PARAF	927091153
GTIN	08023279350913
GMDN	45607

Description

Fingertip pulse oximeter to monitor pulse oxygen saturation, pulse rate, and perfusion index.

- it is light, small in size and easy to take
- accurately measures SpO2 value and pulse rate value
- dual-color display of SpO2, PR, PI, Plethysmogram and pulse bar
- measuring range, SpO2 35 - 99%, Pulse rate 30 - 240 bpm
- innovative 4 directions display
- automatic power On/Off
- more than 30 hours continuous working with 2 AAA alkaline batteries
- audible and visible alarm function
- battery voltage indication, low battery voltage alarm
- prevention against water ingress, drop resistant
- manual in 8 languages: GB, FR, IT, DE, ES, PT, GR, PL.

Size: 66 x 36 x 33 mm