



Product Code	35095
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	1844281
CND	Z1203020408
EAN/UPC	8023279350951
PARAF	980478352
GMDN	45607

Description

Professional oximeter to monitor pulse oxygen saturation and pulse rate through patient's finger (adult, child, neonate).

- wireless Bluetooth technology for transmission to PC Windows (XP and Seven) and Android (need APP). Software for data management
- external pediatric SpO2 probes in 4 sizes for pediatric and neonatal use
- large true colour OLED display
- SpO2, PR, Pulse bar, P.I. and plethysmogram
- accuracy from the first reading in a few seconds
- extra stability against motion during measurement
- audible and visible alarms with high and low alarm limits adjustable by user
- patient oriented display with automatic change of 4 directions
- 9 languages manual (GB, FR, IT, DE, ES, GR, PT, PL, Arabic)

Communication protocols on request
Supplied in nylon bag with batteries

Technical Specifications

- Pulse rate:
 - range: 30-240 bpm
 - accuracy: ± 2 bpm or 2%
- SpO2:
 - range: 70%-100%
 - accuracy: $\pm 3\%$
- Indice di perfusione: 0.2-20%
- Alarm: SpO2: 85-99%
- Puls
 - high 100-240 bpm
 - low 30-60 bpm
- Power: 2xLR03 (AAA) Alkaline batteries 3.0 VDC
- Size: 60 x 33 x 30 mm
- Weight: 35 g