

SpO2 PROBE for PC-3000 and VITAL - paediatric



Product Code	35136
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

CH-REP	Yes
RDM (NSIS)	1729831
CND	Z1203020280
EAN/UPC	8023279351361
GMDN	36087

Description

Accessories - pediatric SpO2 probe spare for PC 3000 Multi-parameter, Vital Signs monitor and UP 7000 (code 35120/21/22/23/24/32 - 35134 - 35145).

For UP 7000 monitor up to s/n J3300RF01475