



Product Code	23514
Unit of sale	1 pc
Minimum order	1
Type	IVD - Self Test
Class	

CH-REP	No
PARAF	950851586
GTIN	00856362005159
GMDN	62537

Description

IHEALTH GLUCO+ BG5s SMART GLUCOSE MONITOR KIT

Pocket sized, ergonomic, wireless glucose monitor with advanced technology and high reliability.

The Wireless Smart Gluco-Monitoring System is a convenient way to measure and record your blood glucose levels from your mobile device.

- view results instantly on the wireless meter
- save up to 500 test results on the wireless meter
- set up reminders and view glucose and insulin trends in multiple ways
- track test strip quantity and expiration dates automatically, no counting or coding required
- works on-line and off-line, with data automatically uploading at each connection.

Kit includes 10 lancets, lancing device, USB charging cable and travel case.

iHealth Gluco-Smart App, available in GB, FR, IT, ES, DE, CZ, SK, HR, TR, IR, CN, Arabic.
It allows to manage all your smart iHealth devices and track relations among parameters.

Technical Specifications

Machine size: 98x35x27.8 mm

Measuring method: Amperometric technology using glucose dehydrogenase (FAD-GDH)

Result range: 20 mg/dL~600 mg/dL (1.1 mmol/L~33.3 mmol/L)

Power source: DC 3.7 V

Blood source: Fresh capillary whole blood

Blood volume: Min. 0.5 micro liter