

EAR THERMOMETER for PC-300, OXY-110 - adult - spare



|               |                |
|---------------|----------------|
| Product Code  | 35170          |
| Unit of sale  | 1 pc           |
| Minimum order | 1              |
| Type          | Medical device |
| Class         | II A           |

|            |               |
|------------|---------------|
| CH-REP     | Yes           |
| RDM (NSIS) | 1792026       |
| CND        | V03010102     |
| EAN/UPC    | 6941900605251 |
| PARAF      | 977817701     |
| GMDN       | 63711         |

Description

Ear infrared thermometer