

# USB CABLE for connection PC-300 - GLUCOSE MONITOR 24108,24110,24111,24114



Product Code	35177
--------------	-------

Unit of sale	1 pc
--------------	------

Minimum order	1
---------------	---

Type	No medical device
------	-------------------

Class	
-------	--

UK-REP	Not required
--------	--------------

CH-REP	Not required
--------	--------------

EAN/UPC	8023279351774
---------	---------------

PARAF	987469424
-------	-----------

## Description

USB cable for Gima glucose monitor (24108/10/11/14)