



Product Code	35189
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	Yes
CH-REP	No

RDM (NSIS)	1967409
CND	Z1203020501
EAN/UPC	6934440101130
UMDNS	16173

Description

AIR BP BLOOD PRESSURE UNIT WITH DISPLAY - Wireless

Ultra-portable smart upper arm blood pressure device.
Includes main unit, cuff and charge cable

- one battery charge for 1000 measurements
- easy to use, robust and durable
- professional measurement
- readings:
 - systolic, diastolic, pulse rate
 - pulse pressure, MAP (Mean Arterial Pressure)
 - irregular heartbeat
- works with smartphone
- real-time Virtual Mercury Column in App during measurement
- real-time Pulse Waveform
- multi-patient management
- unlimited recording
- works with Checkme Pro 35186, 35187
- ViHealth App.

Supplied with cuff, charge cable and multilingual manual (GB, FR, IT, ES, DE).

Technical Specifications

Cuff: Upper arm circumference: 22-42 cm
 Connection: Wireless
 Charge input: Micro USB, DC 5V
 Battery life: 1000 measurements, rechargeable
 Pressure measurement range: 0-300 mmHg \pm 3 mmHg
 Pulse rate: 40 to 200 bpm \pm 2 bpm
 Mobile App for iOS: iOS 9.0 or above, iPhone 4s/iPad 3 or above
 Mobile App for Android: Android 4.4 or above, with Bluetooth 4.0 BLE
 Size - weight: 68x20 mm (main unit) - 30g (main unit)
 Clinical accuracy: meet IEC 80601-2-30