

SpO2 PEDIATRIC PROBE for BCI - 0.9 m cable



Product Code	35220
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II B

CH-REP	No
--------	----

RDM (NSIS)	2621198
CND	C900301
EAN/UPC	6900234100139

Description

OEM - Brand: BCI

N: 5

Sensor type: Pediatric finger

Cable Length: 0.9 m

Extension cable: (code 35300) *This SpO2 sensors need extension cable for connection to device.

All products are supplied with manual in GB, FR, IT, ES.

Compatibility: BCI: 3401, 3304, 3303, 3302, 3301, 3300; Spacelabs Burdick: PRO2; Draeger: Oxipac 2000/2500