

SpO2 NEONATAL "Y" PROBE for CSI - 3.0 m cable



Product Code	35225
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II B
CH-REP	No
EAN/UPC	6900234100535

Description

OEM - Brand: CSI (Criticare)
N: 1
Sensor type: Neonatal Y-type
Cable Length: 3.0 m
All products are supplied with manual in GB, FR, IT, ES.

Compatibility: Alaris (IVAC): VitalCheck 4400; Criticare: 503DX. 503 Spot, 504DX/P. 504DX, series 506 (except 506LN3), 507ELC2, 8100 Poet Plus, Poet (602-3A), Scholar, Scholar III