

SpO2 ADULT PROBE for CSI - 3.0 m cable



Product Code	35227
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II B

CH-REP	No
--------	----

RDM (NSIS)	2621198
CND	C900301
EAN/UPC	6900234100276

Description

OEM - Brand: CSI (Criticare)

N: 1

Sensor type: Adult Finger

Cable Length: 3.0 m

All products are supplied with manual in GB, FR, IT, ES.

Compatibility: Alaris (IVAC): VitalCheck 4400; Criticare: 503DX. 503 Spot, 504DX/P. 504DX, series 506 (except 506LN3), 507ELC2, 8100 Poet Plus, Poet (602-3A), Scholar, Scholar III