

SpO2 ADULT PROBE for DATASCOPE - 0.9 m cable



| | |
|---------------|----------------|
| Product Code | 35230 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | Medical device |
| Class | II B |

| | |
|--------|----|
| CH-REP | No |
|--------|----|

| | |
|------------|---------------|
| RDM (NSIS) | 2621198 |
| CND | C900301 |
| EAN/UPC | 6900234100283 |

Description

OEM - Brand: Datascope

N: 1

Sensor type: Adult Finger

Cable Length: 0.9 m

Extension cable: (code 35301) *This SpO2 sensors need extension cable for connection to device.

All products are supplied with manual in GB, FR, IT, ES.

Compatibility: Accusat, Accutorr 3/4 SAT, Accutorr Plus, Passport (Database SpO2)