

SpO2 ADULT PROBE for DATEX-OHMEDA - 0.9 m cable



Product Code	35232
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II B

CH-REP	No
--------	----

RDM (NSIS)	2621198
CND	C900301
EAN/UPC	6900234100306

Description

OEM - Brand: Datex-Ohmeda

N: 1

Sensor type: Adult Finger

Cable Length: 0.9 m

Extension cable: (code 35302) *This SpO2 sensors need extension cable for connection to device.

All products are supplied with manual in GB, FR, IT, ES.

Compatibility: Datex-Ohmeda: AS/3 ,CS/3, CH/S, CH/2S, CH-RS, SAF-F, CG/2S, Cardiocap I, Cardiocap II, Cardiocap 5, Capnomac Ultima, Light Monitor, M-Estpr, Oscar II, OSP-200, Satlite Trans, Satlile, Trans Plus, System 5, ULTS00