

# SpO2 ADULT PROBE for DATEX-OHMEDA - 3.0 m cable



|               |                |
|---------------|----------------|
| Product Code  | 35234          |
| Unit of sale  | 1 pc           |
| Minimum order | 1              |
| Type          | Medical device |
| Class         | II B           |

|            |               |
|------------|---------------|
| CH-REP     | No            |
| RDM (NSIS) | 2621198       |
| CND        | C900301       |
| EAN/UPC    | 6900234100320 |

## Description

OEM - Brand: Datex-Ohmeda  
N: 1  
Sensor type: Adult Finger  
Cable Length: 3.0 m  
All products are supplied with manual in GB, FR, IT, ES.

Compatibility: Datex-Ohmeda: AS/3 ,CS/3, CH/S, CH/2S, CH-RS, SAF-F, CG/2S, Cardiocap I, Cardiocap II, Cardiocap 5, Capnomac Ultima, Light Monitor, M-Estpr, Oscar II, OSP-200, Satlite Trans, Satlile, Trans Plus, System 5, ULTS00