

SpO2 ADULT "Y" EAR PROBE for GE DATEX-OHMEDA - 3.0 m cable



Product Code	35241
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II B

CH-REP	No
RDM (NSIS)	2621198
CND	C900301
EAN/UPC	6900234100559

Description

OEM - Brand: GE Datex-Ohmeda

N: 1

Sensor type: Adult Y-type ear

Cable Length: 3.0 m

All products are supplied with manual in GB, FR, IT, ES.

Compatibility: GE (Datex-Ohmeda): S/5, i4 Modular Monitors with Eseries modules, Corometrics® 259 with TruSignal Ohmeda Technology