

# SpO2 ADULT PROBE for GE - 3.0 m cable



Product Code	35245
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II B

CH-REP	No
--------	----

RDM (NSIS)	2621198
CND	C900301
EAN/UPC	6900234100368

## Description

OEM - Brand: GE  
N: 1  
Sensor type: Adult Finger  
Cable Length: 3.0 m  
All products are supplied with manual in GB, FR, IT, ES.

Compatibility: GE: Dash 2000/3000/4000 with Nellcor option, Dash 3000/4000/5000 with Ohmeda option, Solar 8000/8000M, i/9500 with Tram Module 451/451N/851/851N, CardioLab®/Mac-Lab® with Tram Module 451/451N Corometrics® 120/2120 series, Eagle 1000/3000/4000 Responder 3000 with Nellcor option