

SpO2 ADULT PROBE for GE TRUSIGNAL - 0.9 m cable



Product Code	35247
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II B

CH-REP	No
--------	----

RDM (NSIS)	2621198
CND	C900301
EAN/UPC	6900234100375

Description

OEM - Brand: GE TruSignal

N: 1

Sensor type: Adult Finger

Cable Length: 0.9 m

Extension cable: (code 35303) *This SpO2 sensors need extension cable for connection to device.

All products are supplied with manual in GB, FR, IT, ES.

Compatibility: TruSignal sensors and interconnect cables are backwards compatible with Datex, ohmeda and GE SpO2 technologies.