

SpO2 ADULT PROBE for MINDRAY - 3.0 m cable



Product Code	35252
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II B

CH-REP	No
--------	----

RDM (NSIS)	2621198
CND	C900301
EAN/UPC	6900234100399

Description

OEM - Brand: Mindray
N: 1
Sensor type: Adult Finger
Cable Length: 3.0 m
All products are supplied with manual in GB, FR, IT, ES.

Compatibility: PM-600, PM-5000, PM-6000, PM-7000, PM-8000, PM-9000, MEC-1000 (Nellcor based) PM-50