

SpO2 ADULT "Y" EAR PROBE for NELLCOR - 3.0 m cable



Product Code	35261
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II B

CH-REP	No
RDM (NSIS)	2621198
CND	C900301
EAN/UPC	6900234100573

Description

OEM - Brand: Nellcor
N: 3
Sensor type: Adult Y-type ear
Cable Length: 3.0 m
Extension cable: (code 35304) ** Available extension cable.
All products are supplied with manual in GB, FR, IT, ES.

Compatibility: Artema: Diascope, Datascope: (Mindray): Accutorr Plus, Expert, Passport, PM-50; Drager: Cicero, Narkomed, Vitalert 1000, GE Medical: Corometrics 118, 120, 129, 188, 510, 556, Critikon 9710, 9720, Plus, Compact, Dinamap Select, MPS, Oxys Shuttle, Pro 300/400, Eagle Series, Dash Series, Hellige Eagle, 1000, 3000, Eagle SMK, SMU, SMV, SOLAR Series, TRAM Modular x50SL Series; Invivo Research: MI2 (Nellcor SpO2), Millennia 4500 Plus 2, Invivo/MDE Escort 100 OPT11A Escort 300A, Escort II OPT30, Prism (Nellcor SpO2); Ivy Biomedical: 405A, 4050, 405P, 405T, Vital-Guard 450 CN, Mennen: Horizon XL (Nellcor SpO2), MR 1300 Mercury 1200; Nellcor: N100 (req. preamp cable), N180, N200 (req. preamp cable), N-20, NPB-75, NPB-290, NPB-295, N-390, N-395, N-3000, N6000; Nihon Kohden: BSM-1102 Life Scope EC, BSM2304 iPro (Life Scope i), BSM-4102, 4112 Procyon LT (Life Scope P), BSM-4104, 4114 Procyon (Life Scope P), BSM-9510 Life Scope M (Nellcor SpO2), OGS2000 PocketCare (Nellcor SpO2); Pace Tech Vitalmax 800 Series (cara negra Nellcor SpO2): Philips 78352C, 78354C, 78834C, MI020A, M1025B, M1205A, M1350B, M1722A/B, M1732A/B, M2475B, M2601A, M3000A